



Motorized Clinical Microscope with LED Illumination



with Motorized Operation

The Ci-E is equipped with Tilt/Telescopic Binocular Tube and Dual Position Low Profile stage. It is ideally suited for applications and sample analysis that require frequent magnification switching.

Motorized Nosepiece

Motorized Swing-Out Condenser

Top lens element automatically swings-in/out depending on the objective lens that's selected, using intelligent linking.

Auto Light Intensity Reproduction

The user-defined light intensity for each objective lens is automatically memorized and replicated when the objective is used again. This eliminates the need for manual re-adjustment of light intensity.

Nosepiece Rotating Buttons



The nosepiece can be easily rotated allowing you to keep your eyes on the specimen. Simply toggle up or down in magnification with a touch of your finger tips. Or with one press of the button, you can quickly switch between your two favorite magnifications (objective lenses). This is useful when frequent change of magnifications is necessary, for example between 10x and 40x objectives.

Eco-Illumination

Image Capture Button

For more information, go to www.microscope.healthcare.nikon.com/ci or call 1-800-52-NIKON







